

LMA | Classic Excel™

The Original LMA™ Airway with Extended Benefits – Better Than Ever

Used in millions of cases by thousands of clinicians worldwide, LMA™ now brings you the LMA Classic Excel™. Including all of the features of the LMA Classic™, this improved design now facilitates intubation and is reusable up to 60 times.

This enhanced design is ideal for unforeseen airway complications where intubation becomes necessary.

LMA™ airways are the undisputed choice of medical professionals in more than 35 countries. Over 2,500 published clinical references and more than 150 million patients attest to the success of the one and only LMA™ airway family. The LMA Classic Excel™ is more cost efficient and better than ever.


- Reusable up to 60 times
- Removable connector for easy fiberoptic assisted intubation
- Facilitates intubation for up to a size 7.5 ETT
- Epiglottic Elevating Bar (EEB) designed to aid intubation
- More cost effective
- Soft silicone cuff reducing the likelihood of throat irritation and stimulation
- Latex free
- Compared to an ET Tube:
 - Facilitates smooth emergence from anesthesia
 - Less invasive
 - Less drug use
 - Minimal hemodynamic response
- Part of ASA difficult airway algorithm



LMA | Classic Excel™



LMA Classic Excel™ with removable connector allows up to a size 7.5 ETT



LMA Classic Excel™ with an ETT inserted

Sizes & Catalog Numbers

Product Description	Patient Size	Mask Size	Catalog #
LMA Classic Excel™ - 1 Unit	Children 30 – 50 kg	3	902030
LMA Classic Excel™ - 1 Unit	Adults 50 – 70 kg	4	902040
LMA Classic Excel™ - 1 Unit	Adults 70 – 100 kg	5	902050
LMA Classic Excel™ removable connectors – Size 3, 4 (10 pack)	N/A	N/A	904900
LMA Classic Excel™ removable connectors – Size 5 (10 pack)	N/A	N/A	904901
Cuff-Deflator (Blue)	N/A	3, 4	99130
Cuff Deflator (Gray)	N/A	5	99150
BiteGard™ (50 single use bite blocks)	N/A	N/A	99300
LMA Classic Excel™ Instructions for Use	N/A	N/A	Available on our website or on request
LMA™ Airway Instructional CD-Rom	N/A	N/A	Available on our website or on request

Selection Guidelines

Mask Size	Patient Size	Maximum Cuff Inflation Volume	Maximum Size ETT
3	Children 30 – 50 kg	Up to 20 mL	7.0 mm Cuffed*
4	Adults 50 -70 kg	Up to 30 mL	7.0 mm Cuffed*
5	Adults 70 – 100 kg	Up to 40 mL	7.5 mm Cuffed*

*LMA Reusable ET Tube is larger than the average ETT – compatible sizes for the LMA Reusable ETT are: LMA Classic Excel™ size 3 – 6.5 mm Cuffed ETT, LMA Classic Excel™ size 4 – 6.5 mm Cuffed ETT, LMA Classic Excel™ size 5 – 7.0 mm ETT Cuffed.

To place an order call 1-800-788-7999 or fax 858-622-4130



Better by Design

4660 La Jolla Village Drive, Suite 900 • San Diego, CA 92122 • 800-788-7999 • Fax 858-622-4130 • www.LMANA.com

Consult the instructions on indications, contraindications, warnings and precautions, or information on which LMA™ airways are best suited for different clinical applications. LMA and LMA Classic Excel are trademarks of The Laryngeal Mask Company Limited. © 2008 LMA North America, Inc. All rights reserved. LMA-561 8/08